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In case of traumatic and degenerative diseases of joints there is a reduction in the quantity of hyaluronic acid in situ, with a consequent reduction of the viscosity of the synovial fluid. This shortage compromises the functionality of joints and causes a series of painful symptoms. The intra-articular injection of hyaluronic acid is able to restore the visco-elastic properties of the synovial fluid, with a consequent relief of pain an evident anti-inflammatory response and an improvement of articular mobility.

The recommended **REGENFLEX®** protocol consists of a two-stage intervention: the first one involves 2 weekly injections of **REGENFLEX® STARTER** 1,6%, whereby the target is to favourably condition the joints thanks to impact, which reduces the symptoms of the specific disease, and maintains a prolonged permanence in situ (in virtue of higher concentration). The following injections of **REGENFLEX®** 0,8% (one weekly session for 3 weeks) will be then carried out in order to guarantee a maintenance, and stabilizing function. **REGENFLEX®** presents optimal tolerability, as it only acts locally (influencing only the joint it is injected to), thus excluding possible systemic actions. Clinical trials have shown how **REGENFLEX®** determines a general improvement of the interested area, with a reduction of pain and swelling, and, functional recovery of the joint. These advantages can last for up to 24 weeks.

#### INSTRUCTIONS FOR USE

##### **REGENFLEX® STARTER**

1 pre-filled disposable syringe (32,0 mg of acido ialuronico sale sodico in 2 ml of buffered physiological solution). Steam-sterilized syringes.

##### **REGENFLEX®**

1 pre-filled disposable syringe (16,0 mg of hyaluronic acid sodium salt in 2 ml of buffered physiological solution). Steam-sterilized syringes.

##### ► Composition

Main ingredient: Hyaluronic acid sodium salt.

Other ingredients: Sodium chloride, Phosphate buffer, water for injections.

##### ► Indications

**REGENFLEX®** e **REGENFLEX® STARTER** are substitutes for synovial fluid, capable of restoring the physiological and rheological properties of compromised joints, in case of pain or reduced mobility caused by degenerative or post-traumatic diseases.

##### ► Dosage and instructions for use

**REGENFLEX® STARTER:** Inject 2 ml of product into the joint once a week for three weeks, subject to the physician's approval and according to health conditions of the patient. Whenever further cycles of treatment are needed, a six-month interval is recommended. Depending on the doctor's opinion and the patient's health conditions, more than one joint may be treated at a time.

**REGENFLEX®:** Inject 2ml of product into the joint once a week for five weeks, subject to the physician's approval and according to health conditions of the patient.

Whenever further cycles of treatment are necessary, a six-month interval is recommended. Depending on the doctor's opinion and the patient's health conditions, more than one joint may be treated at a time.

The volume scale on the syringes is merely indicative.

##### ► Usage

Remove by aspiration any articular effusion before proceeding to the injection of **REGENFLEX®** or **REGENFLEX® STARTER**. Carefully remove the cap from the syringe in order to avoid any contact with the opening. Insert a proper needle (from 18 to 22 G), twisting it on tightly so to ensure that it is airtight and to prevent leakage of the solution.

Inject the product only inside the synovial space, at room temperature and in conditions of complete asepsis.

##### ► Warnings and contraindications

Do not use after the expiry date printed on the box.

Do not use if the box is damaged or open.

The skin where the injection is to be carried out must be sound.

Do not carry out intra-vascular injections.

Do not inject into the knee joint when there are venous or lymphatic stasis of the limb present.

Do not inject outside the joint space, into the synovial tissue or into the joint capsule.

Do not administer in presence of excessive intra-articular effusions.

**REGENFLEX® STARTER** and **REGENFLEX®** must not be injected in case of cutaneous inflammations, infections or pathological diseases in the area to be treated.

After intra-articular injection the patient should be advised to avoid intensive physical activity and to take up normal activities only after a few days.

The content of the pre-filled syringe is sterile.

The external surface of the syringe is not sterile.

Do not sterilize again. The device is disposable.

Store at a temperature lower than 25 °C and away from heat sources. Do not freeze.

Once opened the device must be used immediately and eliminated after use.

Keep out of the reach of children.

##### ► Side effects

Extra-articular injections may cause local side effects.

Where the needle is inserted there can be symptoms like pain, heat, reddening and swelling. In this case it is useful to apply ice to the joint treated. In normal conditions these symptoms disappear in a short time.

The doctor must ensure that the patient informs him of any side effects that may eventually occur after the treatment.

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In rare cases the product can cause allergic reactions.

**REGENFLEX® STARTER** AND **REGENFLEX®** ARE ONLY SOLD UPON A MEDICAL PRESCRIPTION.  
INTRA-ARTICULAR INJECTION MUST ONLY BE CARRIED OUT BY A DOCTOR.

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